



## TR CLASSIC TILT TABLE

---

The TR Classic Tilt Table provides a smooth and effortless lift from horizontal to upright position.

The TR Classic Tilt Table is designed to make testing, rehabilitation, positioning and transferring patients easier and safer. The TR Classic Tilt Table is a logical choice for patients who have paralysis, are bedridden or wheelchair bound who would benefit from passive standing from improving circulation to reduction of lower extremity edema. The TR Classic Tilt Table can be used in the home to assist the caregiver with transfer and treatments needs. The TR Classic Tilt Table operates on a standard 115V or 230 V power source for smooth and silent operation. Standard equipment includes 3 safety belts for comfort and security, footboard, removable wire basket and hand operated remote control. Optional activity tray available.

*... innovations in care*



## Tilt Table Features

The TR Classic Tilt Table is electrically operated for smooth and silent operation. It effortlessly tilts from horizontal to a 90 degree upright position.

A hand control is provided for convenient tilt adjustment and tilt indicator for precise monitoring of positioning.

The wide footplate, adjustable and lockable in three angle positions, meets the varied needs of the patient population.

The table comes standard with a back-up battery for emergency lowering in case of a power failure. Wire storage basket is removable.

The TR Classic Tilt Table easily maneuvers on four 12,5 cm (5") locking casters with straight steering. Wheel fenders protect the lying surface when transporting through narrow doorways.

TR reserves the right to make design and/or improvements without notice.

## MODELS:

Standard 0647500  
Pediatric 0647510

## SPECIFICATIONS:

### Dimensions:

Cushion length-standard: 175 cm (69")  
Cushion length-pediatric: 140 cm (57")  
Cushion width-standard: 65 cm (25-1/2")  
Cushion width-pediatric: 45 cm (17-3/4")  
Height horizontal position: 600 cm (23-5/8")  
Casters: 12,5 cm (5")

### Material:

Frame: powder coated steel  
Cushion: foam filled, vinyl cover

### Power:

Electric: 115V/230V  
Emergency lowering battery backup: 9 volt

### Weight Capacity:

150 kg (330)

### ACCESSORIES:

Activity Tray

CE mark



TRANÅS ROSTFRIA AB  
P.O.Box 116  
SE-573 22 Tranås, Sweden  
Phone +46 (0)140 385060  
Fax +46 (0)140 385061  
[www.trequipment.se](http://www.trequipment.se)

... innovations in care