



FOLDABLE SHOWER TROLLEY






FOLDABLE SHOWER TROLLEY

The Leit foldable shower trolley provides a lightweight, folding changing platform for when space is at a premium in the bathroom. It is ideal for use on holiday or travelling as it can be easily folded and transported. Made from aluminium, the trolley weighs approximately 18.5kg and available in 2 sizes; 1400mm or 1640mm. It can be used in conjunction with mobile hoists (please enquire for compatibility).

It can be used with the Leit shower changer or be supplied with a vinyl folding top section. Two side rails are also provided as standard.




	A	B	C	D	E	F	G
Model no.	Length	Width	Deck Height	Folded Length	Folded Width	Folder Height	Ground Clearance
	1714mm	705mm	852mm	237mm	705mm	926mm	1200mm

Approximate dimensions

User Weight: 120kg, Unit Weight: 18.5kg, Maximum Load: 120kg

We reserve the right to alter specifications for the process of continuous development without notice



standard features

- Two side rails
- Folding base
- Folding top section 1400 and 1670 mm
- Lockable castors



additional options

- Travel bag
- Head/foot rest extender
- Compatible with Leit Shower Changer