



# TR Classic Tilt Table



The TR Classic Tilt Table is designed to make testing, rehabilitation, positioning and transferring patients easier and safer.

It is also a logical choice for patients who benefit from passive standing to improve circulation and reduce lower extremity edema.

It can be used in all environments to assist the caregiver with transfer and treatment needs.

TR Classic Tilt Table  
RRP : \$4,650.00



Special price  
\$2,990.00 +  
Freight

Dont miss it!  
Contact Paul McKay @ AIS Healthcare for  
demonstration or further information

Paul McKay  
National Business Development Manager  
email: paul@aispl.com.au  
m: 0407 077 156

## Specifications

Cushion length - Standard	1750 mm (Stocked item)
Cushion length - Pediatric	1400 mm (Non stocked item)
Cushion width - Standard	650 mm (Stocked item)
Cushion width - Pediatric	450 mm (Non stocked item)
Height (Horizontal position)	600 mm
Working Load Limit	150 kg
Power supply	230 V 50 Hz
Power failure lowering	Battery backup (9 V)
Frame	Powder coated steel
Castors	125 mm Two castors with brakes and two with straight steering
Control	Hand control
Cushion	Foam filled, vinyl cover



Registered at Medical Products Agency, Sweden,  
as a MD Class 1 product, MDD 93/42EEC.



## Features



### Rehabilitation

Logical choice for rehabilitation patients who benefit from passive standing to improve circulation and reduce lower extremity edema



### Tilts effortlessly

Tilts effortlessly from horizontal to a 90 degree upright position and tilt indicator for precise position monitoring



### Easy manoeuvring

Easy manoeuvring thanks to large, locking castors with straight steering



### Safety belts

Three safety belts ensure resident safety. Adjustable positioning of the belts makes tilt table suitable for many different residents



### Wide footplate

Wide footplate, adjustable and lockable in three angle positions for resident comfort and safety



### Battery backup

Built-in battery backup for lowering in case of a power failure



### Operation

Smooth and silent operation thanks to an electrical actuator



### Hand control

Hand control for convenient tilt adjustment



### Activity tray

Activity tray can be attached to use for reading or other activities